



**CAN LOCK
CHAIN TYPE CLAMP
MODEL CLL-010-CS
For Lab Mixers**

Features:

- Rod Block with wing nut
- Block bore for 1" diameter rod
- Chain type clamp
- Insertion rod, 3/8" diameter
- Chain, 15" long
- Chain tension adjustment
- Plastic cover on fork holding the mixing can

This can lock may be used on all Dantco Laboratory type Dispersers or Mixers in Series 501, 502 and 503.

When the mixer stand vertical rod is secured in the rod block, the can lock will hold any mixing container with a diameter of 3" to 7". The chain tension is adjusted to lock the container in place.

The can lock has only two components, Chain Clamp (P/N 539-2432-001) and Rod Block (P/N 539-2433-001).



DANTCO MIXERS CORPORATION
9 OAK STREET
PATERSON, NJ 07501-3013
973-278-8776, FAX 973-278-1538
EMAIL sales@dantco.com, WEB <http://www.dantco.com>